



OMRON

AUTOMATIC BLOOD PRESSURE MONITOR

HEM-7322

Comfortable, Advanced Functionality
Providing accurate and comfortable upper arm blood pressure measurement with enhanced IntelliSense Technology

ENHANCED

INTELLISENSE TECHNOLOGY
Automatically applies the right amount of pressure for fast, accurate and more comfortable measurements.



Cuff Wrapping Guide Lamp



Average of Last 3 Readings



IntelliSense Technology



Blood Pressure Colour Indicator



Body Movement Detection



Blood Pressure Level Indicator



90 Memory with Date and Time



Irregular Heartbeat Detection



AC ADAPTER INCLUDED



PREMIUM

All for Healthcare

TRUSTED WORLDWIDE

OMRON HEM-7322 is a compact, fully-automatic blood pressure monitor that offers users a personalised blood pressure monitoring experience while promising accuracy, comfort and ease of use.

Intelli sense ENHANCED INTELLISENSE TECHNOLOGY

The new algorithm enhances measurement accuracy with a standard deviation of less than 4 mmHg. With a press of a button, the enhanced IntelliSense Technology delivers accurate measurements with comfort and ease of use.

OK CUFF WRAPPING GUIDE LAMP

The "Cuff Wrapping Guide Lamp" is a unique feature that guides user to the right cuff wrap by indicating "OK". This feature ensures an accurate and reliable blood pressure measurement.

- : Cuff has been correctly wrapped.
- : Cuff is incorrectly wrapped. Rewrap cuff.

SYS DIA BLOOD PRESSURE COLOUR INDICATOR

For easy interpretation of blood pressure results. When either the SYS or DIA readings are lower than 135 systolic/85 diastolic, the corresponding LED glows in Green. If the SYS or DIA readings are higher than or equal to 135 systolic/85 diastolic, the corresponding LED glows in Orange.

135 85 BLOOD PRESSURE LEVEL INDICATOR

For easy interpretation of blood pressure results. Recent studies* suggest that systolic or diastolic pressure values above 135 systolic/ 85 diastolic mmHg indicate high blood pressure. This criteria is for home blood pressure measurement.

*Source: American Heart Association (AHA) guidelines for home BP monitoring, Ohasama Study



DETECTS IRREGULAR HEARTBEAT

An irregular heartbeat is a heartbeat rhythm that varies by more than 25% from the average heartbeat rhythm detected while the unit is measuring the systolic and diastolic blood pressure.

If the irregular heartbeat symbol is shown frequently, please inform your doctor.



AVERAGE OF LAST 3 READINGS TAKEN WITHIN 10 MINUTES



DETECTS BODY MOVEMENT ERROR

Prompts user to retake measurement when an error is detected due to movement of the body during measurement.



STORES 90 MEASUREMENTS WITH DATE AND TIME IN MEMORY

SPECIFICATIONS

Model	HEM-7322
Display	LCD Digital Display
Measurement Method	Oscillometric method
Measurement Range	Pressure: 0 to 299 mmHg • Pulse: 40 to 180 beats/min.
Accuracy	Pressure: ±3 mmHg • Pulse: ±5% of display reading
Inflation	Fuzzy-logic controlled by electric pump
Deflation	Automatic pressure release valve
Power Source	4 "AA" batteries 1.5V or AC adapter (INPUT AC100-240V 50/60Hz 0.12A)
Battery Life	Capacity of new alkaline batteries is approximately 1000 measurements
Console Weight	Approximately 380g (without batteries)
Cuff Weight	Approximately 130g
Outer Dimensions	Approximately 124(W)mm x 90(H)mm x 161(L)mm
Cuff Dimensions	Approximately 145mm x 466mm (fits arm circumference 22 to 32cm)
Cuff Material	Nylon, polyester, polyvinyl chloride (latex-free)
Package Contents	Main unit, M-size Soft Cuff, AC adaptor, battery set, storage case, instruction manual
Optionals	Wide Range Soft Cuff HEM-RML31 (fits arm circumference 22-42cm) M-size Soft Cuff HEM-CR24 (fits arm circumference 22-32cm) S-size Soft Cuff HEM-CS24 (fits arm circumference 17-22cm) • AC adapter "S"



*Note: Specifications may be subject to change without prior notice.

OMRON HEALTHCARE Co., Ltd.
53, Kunotsubo, Terado-cho
Muko, Kyoto, 617-0002 JAPAN

Distributed by:

Asia Pacific HQ
OMRON HEALTHCARE SINGAPORE PTE LTD.
438A Alexandra Road, #05-05/08
Alexandra Technopark, Singapore 119967
www.omronhealthcare-ap.com

EC REP
OMRON HEALTHCARE EUROPE B.V.
Scorpius 33, 2132 LR Hoofddorp
THE NETHERLANDS

CE0197