

OMRON

Digital Thermometer

For oral, rectal and underarm temperature measurement

MODEL MC-343F



No Pollutants

10 seconds[#]
Rectal

20 seconds[#]
Oral

- ❖ Professional accuracy
- ❖ **Flexible Tip for added comfort and safety**
- ❖ Memory
- ❖ With Buzzer
- ❖ Water-resistant
- ❖ Digital, no fear of mercury poisoning




Omron digital thermometer offers a safe, accurate and quick temperature reading. It eliminates any worries about broken glass or mercury hazards.

For oral, rectal and underarm temperature measurement

- ❖ Quick actual measurement
 - #20 sec. for oral, 10 sec. for rectal and 25 sec. for underarm measurement
 - *25 seconds is the earliest estimated temperature indication. To complete the measurement, it is required to measure the temperature in the underarm between 2 - 5 minutes (depending on the ambient temperature and condition)
- ❖ Flexible tip provides added comfort and safety
- ❖ Accurate measurement- measures actual temperature in °C and °F
- ❖ Water-resistant for easy cleaning
- ❖ Memory- recalls last reading
- ❖ Battery can be replaced
- ❖ Buzzer- Beeps when measurement is completed


Specifications:

Model	MC-343F
Power supply	1.5V DC, 1 LR41 Alkaline-Magnesium button battery
Power consumption	0.1mW
Battery life	Approx. 2 years or more (3 times /day)
Sensing unit	Thermistor
Measurement method	Actual measurement
Temperature display	3-digit, +°C (°F) display in 0.1 degree increment
Measurement Accuracy	±0.1°C (32.0 to 42.0 °C); ±0.2°F (89.6 to 107.6 °F) when measured at a standard room temperature of 23 °C (73.4 °F) in a test tank held at constant temperature
Measurement Range	32.0 °C to 42.0 °C (89.6 to 107.6 °F)
Weight	12g (with battery)
External Dimensions	19.4mm (w) x 132.5mm (l) x 10.0mm (d)
Memory Function	Recalls previous reading
Auto off function	Approx. 30 minutes after use or 3 minutes when not been used
Battery replacement indication	Battery symbol displayed
Accessories	Protective case, test battery, 5 probe covers, instruction sheet

 **OMRON HEALTHCARE Co.,Ltd**
53, Kunotsubo, Terado-cho,
Muko, Kyoto 617-0002 Japan

Asia Pacific HQ
OMRON HEALTHCARE SINGAPORE PTE LTD
438A Alexandra Road, #05-05/08,
Adexandra Technopark, Singapore 119967
www.omron-healthcare.com.sg

Distributed by:



OMRON HEALTHCARE EUROPE B.V.
Scorpius 33, 2132 LR Hoofddorp
THE NETHERLANDS

CE0197